



Human Superoxide Dismutase 1 ELISA

Catalog # LF-EK0101 (1 kit)

Catalog # LF-EK0102 (4 kits bundle)

Catalog # LF-EK0103 (10 kits bundle)

*Sandwich Enzyme-Linked Immunosorbent Assay for Quantitative Detection of
human SOD 1 (Cu/Zn SOD)*

**For research use only
Not for diagnostic or therapeutic purposes**

Contents

1. Introduction	3
2. Principles of Method	3
3. Intended Use	4
4. Storage and Stability	4
5. Chemical Hazard	4
6. Kit Contents	5
7. Materials Required But Not Provided	6
8. Reagent preparation	6
1) Human SOD1 standard	6
2) Secondary Antibody	6
3) AV-HRP	7
4) Washing buffer	7
5) Sample preparation	7
9. Assay Procedure	8
10. Characteristics	9
1) Typical result	9
2) Sensitivity	10
3) Specificity	10
4) Precision	11
5) Recovery	11
11. Troubleshooting	12
12. Reference	13

1. Introduction

Superoxide dismutase (SOD) is an antioxidant enzyme involved in the defense system against reactive oxygen species (ROS). SOD catalyzes the dismutation reaction of superoxide radical anion (O_2^-) to hydrogen peroxide, which is then catalyzed to innocuous O_2 and H_2O by glutathione peroxidase and catalase. Several classes of SOD have been identified. : Cu/Zn SOD (SOD1) is localized in cytosol, Mn SOD (SOD2) in mitochondria, and ecSOD (SOD3) in extracellular space.

Cu/Zn SOD (SOD1) is a homodimer containing copper and zinc that is found almost exclusively in intracellular cytoplasmic spaces.

The level of SOD1 elevates in response to a wide array of mechanical, chemical and biological messengers such as heat shock, shear stress, and UVB- and X-irradiation. Decreased levels of SOD1 expression can also be triggered by activation of the AP2 transcription factor. A down-regulation of SOD1 has been shown in alveolar type II epithelial cells and lung fibroblasts after exposure to hypoxia.

The genomic organization of SOD1 genes shows striking similarities between species and consists of five exons and four introns. The TATA and CCAAT boxes, as well as several highly conserved GC-rich regions, have been localized in the proximal promoter region. Also, more than 90 different mutations of the SOD1 gene have been associated with amyotrophic lateral sclerosis (ALS), also known as Lou Gehrig's disease.

2. Principles of Method

The design of this assay is based on a sandwich Enzyme-Linked Immunosorbent Assay (ELISA). The microtiter plate provided in this kit has been pre-coated with a monoclonal antibody specific to human SOD1. Samples are pipetted into these wells. Nonbound SOD1 and other components of the sample should be removed by washing, then biotin-conjugated monoclonal antibody specific to SOD1 added. In order to quantitatively determine the amount of SOD1 present in the sample, Avidin conjugated to Horseradish Peroxidase (HRP) should be added to each microplate well. The final step, a TMB-substrate solution added to each well. Finally, a sulfuric acid solution is added and the resulting yellow colored product is measured at 450nm. Since the increases in absorbency is directly proportional to the amount of captured SOD1.

3. Intended Use

The AbFrontier human Superoxide Dismutase-1 (human SOD1) ELISA kit is to be used for the in vitro quantitative determination of human SOD1 in human serum, human plasma, cell lysate and buffered solution. The assay will recognize both native and recombinant human SOD1.

This kit has been configured for research use only and is not to be used in diagnostic procedures.

4. Storage and Stability

All kit components of this kit are stable at 2 to 8°C. Any unused reconstituted standard should be discarded or frozen at -70°C. Standard can be frozen and thawed one time only without loss of immunoreactivity.

5. Chemical Hazard

- Stop solution: This reagent is an irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothing. Wear suitable protective clothing, gloves and eye protection. In the event of contact with eyes or skin, wash immediately with plenty of water.
- Standard protein and 2nd Antibody containing Sodium Azide as a preservative.

6. Kit Contents

Contents	Number	Volume
96 Well Plate	1 (in aluminum foil bag with desiccant)	
Incubation Buffer	1	30 ml
Washing Buffer	2	(20×) 25 ml
Standard Protein	1 Glass vial (lyophilized)	
Standard/Sample Dilution Buffer	1	25 ml
Secondary Antibody	1 Glass vial (lyophilized)	
AV-HRP	1	150 µl
Secondary Antibody/AV-HRP Dilution Buffer	1	25 ml
Substrate (TMB)	1	15 ml
Stop Solution	1	15 ml
Protocol booklet	1	
Plate sealers	2	

① 96 Well Plate

: Human SOD1 microtiter plate, one plate of 96 wells (8 well strips×12).

A plate using break-apart strips coated with a mouse monoclonal antibody specific to human SOD1.

② Standard Protein

: Recombinant human SOD1.

③ Secondary Antibody

: Biotinylated anti human SOD1 antibody.

④ AV-HRP

: Avidin linked Horseradish Peroxidase (HRP, enzyme)

⑤ Substrate (Stabilized chromogen)

: Tetramethylbenzidine (TMB) solution

⑥ Stop Solution

: 1 N solution of sulfuric acid (H₂SO₄)

⑦ Plate sealer

: Adhesive sheet

- Do not mix or interchange different reagents from various kit lots.

7. Materials Required But Not Provided

- ① Microtiter plate reader capable of measurement at or near 450nm.
- ② Calibrated, adjustable precision pipettes, preferably with disposable plastic tips (A manifold multi-channel pipette is desirable for large assays.)
- ③ Distilled or deionized water
- ④ Data analysis and graphing software
- ⑤ Vortex mixer
- ⑥ Polypropylene tubes for diluting and aliquoting standard
- ⑦ Absorbent paper towels
- ⑧ Calibrated beakers and graduated cylinders of various sizes

8. Reagent Preparation

1) Human SOD1 standard

1. Reconstitute the lyophilized Human SOD1 standard by adding 1 ml of *Standard/Sample Dilution Buffer* to make the 10 ng/ml standard stock solution. Allow solution to sit at RT for 5 minutes, then gently vortex to mix completely. Use within one hour of reconstituting.

2. Prepare 1 ml of 800 pg/ml top standard by adding 80 μ l of the above stock solution in 920 μ l of *Standard/Sample Dilution Buffer*. Perform 2-fold serial dilutions of the top standards to make the standard curve within the range of this assay (12.5 pg/ml ~ 800 pg/ml) as below. *Standard/Sample Dilution Buffer* serves as the zero standard(0 pg/ml).

Standard	Add	Into
800.0 pg/ml	80.0 μ l of the std.(10 ng/ml)	920.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
400.0 pg/ml	500 μ l of the std.(800 pg/ml)	500.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
200.0 pg/ml	500 μ l of the std.(400 pg/ml)	500.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
100.0 pg/ml	500 μ l of the std.(200 pg/ml)	500.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
50.0 pg/ml	500 μ l of the std.(100 pg/ml)	500.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
25.0 pg/ml	500 μ l of the std.(50 pg/ml)	500.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
12.5 pg/ml	500 μ l of the std.(25 pg/ml)	500.0 μ l of the <i>Standard/Sample Dilution Buffer</i>
0 pg/ml	1 ml of the <i>Standard/Sample Dilution Buffer</i>	

2) Secondary Antibody

100X secondary antibody solution can be made by adding 150 μ l secondary antibody/AV-HRP dilution buffer in the vial.

1. Equilibrate to room temperature, mix gently.

2. Mix 20 μ l *Secondary Antibody concentrated solution (100X)* + 1.98 ml *Secondary Antibody/AV-HRP dilution buffer*. (Sufficient for two 8-well strip, prepare more if necessary) Label as “Working Secondary antibody Solution”.
3. Return the unused *Secondary Antibody concentrated solution* to the refrigerator.

3) AV-HRP

1. Equilibrate to room temperature, mix gently.
2. Mix 20 μ l *AV-HRP concentrated solution (100X)* + 1.98 ml *Secondary Antibody/AV-HRP dilution buffer*. (Sufficient for two 8-well strip, prepare more if needed)
Label as “Working AV-HRP Solution”.
3. Return the unused *AV-HRP concentrated solution* to the refrigerator.

4) Washing buffer

1. Equilibrate to room temperature, mix to re-dissolve any precipitated salt.
2. Mix 0.5 volume *Wash buffer concentrate solution (20X)* + 9.5 volumes of deionized water. Label as “Working Washing Solution”.
3. Store both the concentrated and the Working Washing Solution in the refrigerator.

*** Directions for washing**

1. Fill the wells with 300 μ l of “Working Washing Buffer”.
Let soak for 1 to 3 minutes and then all residual wash-liquid must be drained from the wells by aspiration (taking care not to scratch the inside of the well) or decantation, followed by forceful tapping of the plate on absorbent paper. Never insert absorbent paper directly into the wells.
If using an automated washer, the operating instructions for washing equipment should be carefully followed.
2. Incomplete washing will adversely affects the assay and renders false results.
3. It is recommended to use laboratory tape to hold plate strips to the plate frame while performing the plate washing to avoid strips coming free of the frame.

5) Sample preparation

Blood should be collected by veinpuncture. For plasma samples, blood may be drawn into tubes containing sodium citrate or heparin, EDTA. The serum or plasma should be separated from the coagulated or packed cells by centrifugation. Specimens may be shipped at room temperature and then stored refrigerated at 2-8°C if testing is to take place within one week after collection. If testing is to take place later than one week, specimens should be stored at -20°C. Avoid repeated freeze/thawing.

9. Assay Procedure

- Allow all reagents to reach room temperature before use. Gently mix all liquid reagents prior to use.
 - All standards, controls and samples should be run in duplicate for confirmation of reproducibility.
 - A standard curve must be run with each assay.
 - If particulate matter is present in the analyte, centrifuge or filter prior to analysis.
 - Maintain a consistent order of components and reagents addition from well to well. This ensures equal incubation times for all wells.
- 1) Determine the number of 16-well strips needed for assay. Insert these in the flame(s) for current use (Re-bag extra strips and frame. Refrigerate for further use).
 - 2) Add 300 μ l of *Incubation buffer* to all wells and incubate the plate for 5 minutes at room temperature.
 - 3) Thoroughly aspirate or decant the solution from the wells. Wash wells 2 times (See “Directions for washing”).
 - 4) For the standard curve, add 100 μ l of the standard to the appropriate microtiter wells. Add 100 μ l of the *Standard/Sample Dilution Buffer* to zero wells.
 - 5) Serum and plasma require at least 20 fold dilution in the *Standard/Sample Dilution Buffer*. And add 100 μ l of samples to each wells.
 - 6) Cover the plate with the plate cover and incubate for 2 hours at room temperature.
 - 7) Thoroughly aspirate or decant the solution from the wells. Wash the wells 3 times (See “Directions for washing”).
 - 8) Pipette 100 μ l of “Working Secondary Antibody Solution” into each well.
 - 9) Cover the plate with the plate cover and incubate for 1 hour at room temperature.
 - 10) Thoroughly aspirate or decant the solution from the wells. Wash the wells 3 times (See “Directions for washing”).
 - 11) Add 100 μ l “Working AV-HRP Solution” to each well.
 - 12) Cover the plate with the plate cover and incubate for 30 minutes at room temperature.
 - 13) Thoroughly aspirate or decant the solution from the wells. Wash the wells 3 times (See “Directions for washing”).
 - 14) Add 100 μ l of *Substrate* to each well. The liquid in the wells should begin to turn blue.
 - 15) Incubate the plate at room temperature.
 - Do not cover the plate with aluminum foil, or color may develop.
The incubation time for chromogen substrate is often determined by the microtiter plate reader used. O.D. values should be monitored and the substrate reaction stopped before

O.D. of the positive wells exceeds the limits of the instrument. O.D. values at 450nm can only be read after the Stop Solution has been added to each well.

- Because the *Substrate* is light sensitive, avoid the remained *Substrate* solution prolonged exposure to light.
 - Typically, reaction is stopped 5~10 minutes after treatment of Substrate, but this time can be adjusted as the user desires.
- 16) Add 100 μ l of *Stop Solution* to each well. The solution in the wells should change from blue to yellow.
 - 17) Read the absorbance of each well at 450 nm. Read the plate within 20 minutes of adding the *Stop Solution*.
 - 18) Plot on graph paper the absorbance of the standard against the standard concentration (Optimally, the background absorbance can be subtracted from all data points, including standards, unknowns and controls, prior to plotting.). Draw a smooth curve through these points to construct the standard curve.
 - 19) Read the human SOD1 concentrations for the unknown samples and controls from the standard curve plotted in step 18. Multiply value(s) obtained for the unknown sample by the dilution factor (Samples producing signals greater than that of the highest standard should be further diluted in the *Standard/Sample Dilution Buffer*).

10.Characteristics

1) Typical result


The standard curve below is for illustration only and **should not be used** to calculate results in your assay.

A standard curve must be run with each assay.

Standard Human SOD1 (pg/ml)	Optical Density (at 450nm)
0	0.093
12.5	0.102
25.0	0.126
50.0	0.179
100.0	0.289
200.0	0.502
400.0	0.880
800.0	1.754

< Limitations >

- Do not extrapolate the standard curve beyond the 800 pg/ml standard point.
- Other buffers and matrices have not been investigated.
- The rate of degradation of native human SOD1 in various matrices has not been investigated.



2) Sensitivity

The minimal detectable dose of human SOD1 was calculated to be 12.5 pg/ml, by subtracting two standard deviations from the mean of 10 zero standard replicates (ELISA buffer, S0) and intersecting this value with the standard curve obtained in the same calculation.

3) Specificity

The following substances were tested and found to have no cross-reactivity: human SOD2, SOD3, SOD4, rat SOD1 and mouse SOD1.

4) Precision

① Within-Run (Intra-Assay)

(n=5)

Mean (pg/ml)	SD (pg/ml)	CV (%)
204.16	4.23	2.1
400.74	30.84	7.7
798.24	51.20	6.4

② Between-Run (Inter-Assay)

(n=8)

Mean (pg/ml)	SD (pg/ml)	CV (%)
82.76	7.91	9.6
177.79	12.35	6.9
374.29	18.92	5.1

5) Recovery

Recovery on addition is 99.4~99.5 % (mean 99.5 %)

Recovery on dilution is 98.0~106.1 % (mean 100.7 %)

11. Troubleshooting

Problem	Possible Cause	Solution
High signal and background in all wells	• Insufficient washing	• Increase number of washes • Increase time of soaking between in wash
	• Too much AV-HRP	• Check dilution, titration
	• Incubation time too long	• Reduce incubation time
	• Development time too long	• Decrease the incubation time before the stop solution is added
No signal	• Reagent added in incorrect order, or incorrectly prepared	• Review protocol
	• Standard has gone bad (If there is a signal in the sample wells)	• Check the condition of stored standard
	• Assay was conducted from an incorrect starting point	• Reagents allows to come to 20~30 °C before performing assay
Too much signal – whole plate turned uniformly blue	• Insufficient washing – unbound AV-HRP remaining	• Increase number of washes carefully
	• Too much AV-HRP	• Check dilution
	• Plate sealer or reservoir reused, resulting in presence of residual AV-HRP	• Use fresh plate sealer and reagent reservoir for each step
Standard curve achieved but poor discrimination between point	• Plate not developed long Enough	• Increase substrate solution incubation time
	• Improper calculation of standard curve dilution	• Check dilution, make new standard curve
No signal when a signal is expected, but standard curve looks fine	• Sample matrix is masking detection	• More diluted sample recommended
Samples are reading too high, but standard curve is fine	• Samples contain protein levels above assay range	• Dilute samples and run again
Edge effect	• Uneven temperature around work surface	• Avoid incubating plate in areas where environmental conditions vary • Use plate sealer

12. Reference

- 1) Ruby Leah B. et al. *The copper Chaperone CCS Directly Interacts with Copper/Zinc Superoxide Dismutase* (1998) *The Journal of Biological Chemistry* 273(37):23625-23628
- 2) Steen V. Petersen et al. *The dual nature of human extracellular superoxide dismutase: One sequence and two structures* (2003) *PNAS* 100(24):13875-13880
- 3) Tim D. Oury et al. *Human extracellular superoxide dismutase is a tetramer composed of two disulphide-linked dimmers: a simplified, high-yield purification of extracellular superoxide dismutase* (1996) *Biochem .J* .317:51-57
- 4) Tohru Fukai et al. *Extracellular superoxide dismutase and cardiovascular disease* (2002) *Cardiovascular Research* 55:239-249
- 5) Yoshiaki Furukawa et al. *Oxygen-induced maturation of SOD1: a key role for disulfide formation by the copper chaperone CCS* (2004) *The EMBO Journal* 23:2872-2881
- 6) Richard W. Strange et al. *The Structure of Holo and Metal-deficient Wild-type Human Cu, Zn Superoxide Dismutase and its Relevance to Familial Amyotrophic Lateral Sclerosis* (2003) *J. Mol. Biol.* 328:877-891
- 7) Stefan I. Liochev et al. *Cross-compartment protection by SOD1* (2005) *Free Radical Biology & Medicine* 38:146-147
- 8) Richard A. Weisiger et al. *Superoxide Dismutase* (1973) *The Journal of Biol. Chem.* 248(10):3582-3592

◆ Ordering Information

For orders, please contact :

Young In Frontier Co., Ltd.

Tel : +82-2-2140-3300

Fax: +82-2-2140-3310

E-mail: orders@younginfrontier.com

Address: 11F, Byucksan Digital Valley 5th, Gasan-dong 60-73, Geumcheon-gu, Seoul,
Korea (153-801)

Website: <http://www.abfrontier.com>

Or, your local distributor.

For technical advice, please contact:

E-mail : orders@younginfrontier.com

Website : <http://www.abfrontier.com>